

reference.

	Results (n)				Test performance (%)				
	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV	Accuracy
					[95% CI]	[95% CI]	[95% CI]	[95% CI]	[95% CI]
Overall	580	105	308	21	97 (580/601) [95-98]	25 (105/413) [22-30]	65 (580/888) [62-68]	83 (105/126) [76-89]	68 (685/1014) [65-70]
Age									
< 60 years	362	81	225	15	96 (362/377) [94-98]	27 (81/306) [22-32]	62 (362/587) [58-66]	84 (81/96) [76-90]	65 (443/683) [61-68]
≥ 60 years	218	24	83	6	97 (218/224) [94-99]	22 (24/107) [16-31]	72 (218/301) [67-77]	80 (24/30) [63-91]	73 (242/331) [68-78]
Sex									
Male	272	35	148	12	96 (272/284) [93-98]	19 (35/183) [14-25]	65 (272/420) [60-69]	75 (35/47) [61-85]	66 (307/467) [61-70]
Female	308	70	160	9	97 (308/317) [95-99]	30 (70/230) [25-37]	66 (308/468) [61-70]	89 (70/79) [80-94]	69 (378/547) [65-73]

TP= true positive, TN=true negative, FP=false positive, FN=false negative, PPV= positive